

● **Measurement by Ion Selective Electrodes**

Adrian Vazquez

May 8, 2012

● The world leader in serving science

Why Use Ion Selective Electrodes?

- **Responsive over a wide concentration range**
- **Not affected by color or turbidity of sample**
- **Rugged and durable**
- **Rapid response time**
- **Real time measurements**
- **Low cost to purchase and operate**
- **Easy to use**

Why Use ISE's?

- There are many types of ISE applications:
 - *Ammonia in water or wastewater*
 - *Residual chlorine in water or wastewater*
 - *Chloride in water or wastewater*
 - *Cyanide in water or wastewater*
 - *Fluoride in drinking water*
 - *Nitrate in drinking water or wastewater*

3

ThermoFisher
SCIENTIFIC

Why Use ISE's?

- EPA approved methods
 - *Acidity*
 - *Alkalinity*
 - *Ammonia*
 - *Bromide*
 - *Chloride*
 - *Residual Chlorine*
 - *Cyanide*
 - *Fluoride*
 - *Total Kjeldahl Nitrogen (TKN)*
 - *Nitrate*
 - *Dissolved Oxygen/BOD*
 - *pH*
 - *Sulfide*

4

ThermoFisher
SCIENTIFIC

What Are ISE's?

- Electrodes are devices which detect species in solutions
- Electrodes consist of a sensing membrane in a rugged, inert body




5

ThermoFisher
SCIENTIFIC

How Do ISE's Work?

- If two solutions are separated by an ion-permeable membrane, they will equilibrate:




6

ThermoFisher
SCIENTIFIC

How Do ISE's Work?

- If the membrane is permeable to only one species, a charge quickly develops which opposes further movement
- The charge is proportional to the difference in concentration on the two sides
- The total number of ions that diffuse is very small



7

ThermoFisher
SCIENTIFIC

How Do ISE's Work?


- The reference electrode completes the circuit to the sensing electrode (ISE)
- Reference electrodes have a small leak to establish contact with the sample
- The reference solution (usually KCl) in contact with the reference keeps the reference potential constant

8

ThermoFisher
SCIENTIFIC

The Electrode Measurement System

- The meter measures the voltage potential (mV) difference between the sensing electrode and the reference electrode.
- Sensing electrode potential changes with sample concentration.
- Reference electrode potential does not change to compare sensing potential to.



9

ThermoFisher
SCIENTIFIC

ISE Meters

- ISE meters report concentrations
 - *No manual calibration curves are required*
 - *The potentials of known standards are stored in memory to produce a calibration curve*
 - *The potential of the sample is then compared to the calibration curve to determine the concentration*
- ISE meters generate sophisticated curves which are held in the meter's memory
 - *Run standards*
 - *Run unknowns*
 - *Read results*



10

ThermoFisher
SCIENTIFIC

Types Of Sensing Electrodes

- Glass Membrane
- Solid State
- Liquid Membrane
- Gas Sensing

11

ThermoFisher
SCIENTIFIC

Glass Membrane Electrodes

- pH
- Sodium



12

ThermoFisher
SCIENTIFIC

Glass Membrane Electrodes

- **Example: Sodium**
 - *0.02 detection limit*
 - *Available in combination or half-cell versions*
 - *Ross or Ag/AgCl references*



13

ThermoFisher
SCIENTIFIC

Types Of Sensing Electrodes

- **Solid State**



14

ThermoFisher
SCIENTIFIC

Solid State Electrodes


- Bromide
- Cadmium
- Chloride
- Chlorine
- Copper
- Cyanide
- Fluoride
- Iodide
- Lead
- pH
- Silver
- Thiocyanate

15

ThermoFisher
SCIENTIFIC

Solid State Electrodes

- Sensing element is a solid crystalline material
- Store combination electrodes in dilute standard
- Store half-cell electrodes dry
- Polish solid state electrodes to rejuvenate




16

ThermoFisher
SCIENTIFIC

Solid State Electrodes

- **Example: Chlorine**
 - *EPA approved for drinking water or wastewater*
 - *0.01 detection limit*
 - *Combination electrode*



17

ThermoFisher
SCIENTIFIC

Types Of Sensing Electrodes

- **Liquid Membrane**



18

ThermoFisher
SCIENTIFIC

Liquid Membrane Electrodes


- Ammonium
- Calcium
- Chloride
- Fluoroborate
- Nitrate
- Perchlorate
- pH
- Potassium
- Surfactant

19

ThermoFisher
SCIENTIFIC

Liquid Membrane Electrodes

- Sensing membrane is an ion carrier dissolved in a soft plastic
- Store electrodes in dilute standard for short-term
- Store module dry in vial for long-term



20

ThermoFisher
SCIENTIFIC

Liquid Membrane Electrodes

- **Example: Nitrate**
 - *EPA approved for drinking water*
 - *0.1 ppm detection limit*
 - *Available in combination or half-cell versions*



21

ThermoFisher
SCIENTIFIC

Gas Sensing Electrodes

- **Ammonia**
- **Carbon dioxide**
- **Nitrogen dioxide**
- **Oxygen**




22

ThermoFisher
SCIENTIFIC

Gas Sensing Electrodes: Ammonia

- Gas sensing electrodes work by measuring the pH change caused by diffusion of the gas through a hydrophobic but porous membrane



23

ThermoFisher
SCIENTIFIC

Ammonia ISE Measurement

- Measure from 0.01 ppm to 17,000 ppm ammonia
 - Also use for ammonium and organic nitrogen
- Must adjust sample pH
- pH electrode – Measures change in pH across the membrane
- EPA approved method for waste water
- Must pay attention to manufacturer's instructions

24

ThermoFisher
SCIENTIFIC


Ammonia Electrode Setup



Glass sensing surface protruded relatively to the outer body
(view without gas membrane)



Membrane on outer body without the inner sensing stem
(retaining cap not shown)



Membrane on outer body with the inner sensing stem protruded
(retaining cap not shown)

25

ThermoFisher
SCIENTIFIC

Ammonia Tips

- Tilt the electrode in sample or standard
- Rinse probe by immersion
 - More effective than squirting with rinse bottle
 - Does not disturb layer of inner fill solution between membrane and inner stem
- Conditioning electrode.
 - Exposure to ammonia prior to calibration/use greatly improves performance
 - 15 minutes in 1 mg/L standard (with ISA) greatly improves performance

26

ThermoFisher
SCIENTIFIC

Ammonia Tips

- **Electrode Assembly** – Get rid of air bubbles by tapping on electrode body. Pull gently on cable or move sensor up and down
- **Rinsing Electrode** – Rinse generously between measurements, blot dry. Immerse rather than squirt.
- **Position Electrode at Slant** – This will make easier to view and dislodge bubbles
- **Inspect Membrane Surface** – Remove any bubbles from surface of membrane by wiggling or tapping electrode
- **Electrode Storage** – If storing the electrode overnight or up to 1 week, store in the electrode fill solution
- **Stirring** – Stir vigorously but avoid creating bubbles

27

ThermoFisher
SCIENTIFIC

Classic Ammonia ISE Assembly

- Soak inner body in fill solution for 2 hrs.
- Replace the membrane every 2-4 weeks
- Use tweezers to handle membrane
- Stretch membrane until taut
- Use 2.0 - 2.5 mls of fill solution
- Shake electrode down after assembly
- Pull on the electrode cable to allow fresh supply of fill solution to the membrane but re-calibrate after this is done



28

ThermoFisher
SCIENTIFIC

Classic Ammonia ISE Storage

- 9512BNWP Only
- Store in 10-100 ppm ammonia standard
- Between measurements store in 10 ppm standard with NaOH (ISA)
- Condition electrode in pH 4 buffer for several minutes before starting low level calibration

29

ThermoFisher
SCIENTIFIC

High Performance Ammonia ISE Hints


- Membranes and fill solution are NOT interchangeable between old and new electrodes!
- Do not stretch new membrane
 - Just smooth over and pull tight
 - Can use outer bodies with assembled membranes
- Fill electrode to just below fill level line
- Condition assembled electrode overnight in electrode fill solution
- Store the electrode in electrode fill solution or storage solution (951213)
- Between measurements use 1 ppm standard with 1 ml alkaline reagent
- Use Low Level Ammonia ISA (951210). Do not use standard ISA (951211) for any calibrations or measurements below 10 ppm
- Or use 1 ml of alkaline reagent (951011) per 100 mls sample or standard
- Do not reuse fill solution

30

ThermoFisher
SCIENTIFIC

Ammonia Trouble Shooting Tips

- Make sure you have fresh standards and are using the correct amount of ISA
- Verify the electrode been conditioned
- Check slope and drift of electrode
- If you have a poor slope or drift, change the membrane and fill solution. Recheck slope and drift
- Check electrode inner body. Gently handle inner body
- Stir vigorously but without creating bubbles



31

ThermoFisher
SCIENTIFIC

How Are ISE's Used?

- **Direct measurement**
- **Incremental methods**
- **Titration**
- **Gran methods**

32

ThermoFisher
SCIENTIFIC

Direct Measurement

- Preferred method in most cases:
 - *Many samples with similar backgrounds*
 - *High volume of samples*
 - *Wide range of concentrations*
 - *Easy*
- Calibrate by comparison with known standards
- Read by preparing calibration curve or using ISE meter
- Precision is +/- 2%

33

ThermoFisher
SCIENTIFIC

Direct Measurement


- Calibrate every 2 hours
- Always calibrate with standards that bracket expected concentration range
- Always use at least two standards that are ten fold apart in concentration
- Slope range for monovalent ions: 54-60 mV
- Slope range for divalent ions: 26-30 mV

34

ThermoFisher
SCIENTIFIC

Direct Measurement

- Two-point calibration for linear portion of curve
- Low-level measurements require non-linear multi-point calibration or blank correction
- *Use a meter which calculates the blank from the calibration curve*
 - No separate blank has to be run
 - Equivalent to drawing a smooth curve through the lowest three calibration points and extrapolating to zero concentration





35

ThermoFisher
SCIENTIFIC

How Meters Calibrate Slope

- Slope Segmented
 - Point-to-point calibration
- Auto-Blank for Low Levels
 - Smooth curve through points



36

ThermoFisher
SCIENTIFIC

Incremental Methods

- **Known addition**
- **Analate addition**
- **Known subtraction**
- **Analate subtraction**

37

ThermoFisher
SCIENTIFIC

Known Addition

- **A technique which adds a known amount of the species being measured into the sample**
- **This incremental method eliminates the need for a separate calibration curve for each sample**
- **Single, double, or multiple known additions possible**
- **2-4 times more accurate than direct read**

38

ThermoFisher
SCIENTIFIC

Known Addition

- **Known addition is preferred when:**
 - *Number of samples is small*
 - *Backgrounds vary and cannot be fixed*
 - *Working at very low levels*
 - *Excess complexing agent is present*

39

ThermoFisher
SCIENTIFIC

Practical Considerations


- **Method interferences**
- **Electrode interferences**
- **Temperature effects**

40

ThermoFisher
SCIENTIFIC

Method Interferences

- Many method interferences can be overcome by using Ionic Strength Adjusters
- ISA's are added to samples and standards to maintain constant background
- ISA's minimize ionic strength differences
- ISA's can complex interferences
 - Example Fluoride may be complexed by aluminum
- ISA's can adjust pH to proper range



41

ThermoFisher
SCIENTIFIC

Temperature Effects

- A change in temperature will cause electrode response to shift and change slope
- On average, a 1 °C change in temperature gives rise to a 2% error for monovalent ISE's
- On average, a 1 °C change in temperature gives rise to a 4% error for divalent ISE's

Temp (°C)	Slope (mV)
0	54.2
10	56.2
20	58.2
25	59.2
30	60.2
35	61.2
40	62.1
50	64.1

42

ThermoFisher
SCIENTIFIC

Temperature Effects

- Temperature compensation is possible only if the isopotential point of the electrode is used to adjust the calibration curve
- For most situations make sure standards and samples are at the same temperature

43

ThermoFisher
SCIENTIFIC

Troubleshooting Sequence

- Meter
- Standards
- Reference electrode
- Sensing electrode
- Sample
- Technique



44

ThermoFisher
SCIENTIFIC

Troubleshooting pH Problems

- Troubleshooting ISE Meters
 - Use meter shorting strap
 - Reading should be 0 mV +/- 0.2 mV
 - Use meter self-test procedure
- Troubleshooting Standards
 - Use Fresh Standards for calibration
 - Verify expiration date
 - Verify when standard was opened
 - Stir standards during calibration



45

ThermoFisher
SCIENTIFIC

Making Accurate Measurements

- **Maintain a constant reference potential**
- **Adjust ionic strength**
- **Remove method interferences**
- **Remove electrode interferences**
- **Operate at a constant temperature**
- **Stir standards and samples gently**
- **Bracket sample with standards**

46

ThermoFisher
SCIENTIFIC

Thermo Scientific

- Contact us for any technical questions
- Technical Service: (800) 225-1480
- Web site: www.thermoscientific.com/water

47

ThermoFisher
SCIENTIFIC

Q&A



48

ThermoFisher
SCIENTIFIC